



Digital overhead stirrer up to 20lt - AT-MD 20

Summary

Digital overhead stirrer AT-MD 20 is an instrument stirring mechanically a liquid with maximum viscosity of 50.000 mPas up to 20lt.

Description

It is suitable for routine and continuous works

Digital timer

LCD display

Overload protection and over-heat protection for motor

Adjustable spindle to fix rods.

Ros, support, stirring shaft and clamp are not included with instrument. Choose them among accessories according to application.





Features

Item	605.0150.20
Power	V/HZ 230/50-60
Weight	Kg 5.0
Dimension	mm 90x145x180
Power supply	W 83/120 120
Protection class	IP 43 - 40
Display	Digital
Max stirring capacity	Lt 20
Viscosity	mPas 50000
Torque tension	Ncm 145
Rod diameter	mm 10
Variable speed	Rpm 200-3000