



Stirrer without heating - F20

Summary

Electromagnetic stirrer F20 is a model without heating, suitable for continuous applications

It is a friendly product suggested for microbiological and biochemical fields.

Description

Metal structure painted with anti-acid epoxy powder coating gets the instrument strong and robust, the plate is in stainless steel

Direct current motor is useful to work up to 2 liters





Features

Item	601.0101.20
Power	V/HZ 230/50-60
Weight	Kg 0.7
Dimension	mm 110x115x60
Protection class	IP 41
Stainless shelf	mm 110 x 110
Max stirring capacity	lt 2
Variable speed	Rpm 100-1100
Hole for rod	No
Motor Power	W 10
Stirring bar suggested \varnothing	mm 6x20