



MM - 100 ml

Summary

MM heating mantles were developed by FALC on the base of long years of experience;

Description

MM heating mantles were developed by FALC on the base of long years of experience;

- steel housings, oven coated with acidproof epoxy paint;
- semi – rigid structure that fits matrass shape drastically reducing risks of matrass breaking;
- thermal insulation guaranteed by ceramicfiber ;
- temperature control through a graduated thermoregulator.

All models are provided with ON /OFF switch with relevant green light.



Features

Item	625.0500.01
------	-------------



Temperature	+350
External dimension (LxWxH)	mm 185 x 145
Power	V/HZ 230/50-60
Weight	Kg 1.5
Protection class	IP 30
Flask	ml 100
Precision on the plate	± 5
Heating res. capacity	W 130
Palloni \varnothing mm	64