



---

## Digital water bath with squared tank 12 lt - WB-MD15

### Summary

\*CAPACITY 5 - 15 - 24 lt \* External case in ABS CAPACITY 50 lt \* External case in steel painted with epoxy paint ALL CAPACITY

### Description

\*CAPACITY 5 - 15 - 24 lt

\* External case in ABS CAPACITY 50 lt

\* External case in steel painted with epoxy paint ALL CAPACITY

\* Stainless steel tank (made by mould)

\* General switch

\* Digital thermoregulator

\* Safety thermostat

\* Good insulation between external case and tank



## Features

Item	610.1006.15
Temperature	+5 amb. to +120
Precision	°C ±0.5
Heating power	W 1400
Internal dimension (LxWxH)	mm 255x290x150
External dimension (LxWxH)	mm 445x385x240
Power	V/Hz 230/50-60
Weight	Kg 6.2
Security class	3.1
Protection class	IP 42
Capacity	lt 12
Unload + level	No

