



Digital water bath with squared tank 18 lt - WB-MD24

Summary

*CAPACITY 5 - 15 - 24 lt * External case in ABS CAPACITY 50 lt * External case in steel painted with epoxy paint ALL CAPACITY

Description

*CAPACITY 5 - 15 - 24 lt

* External case in ABS CAPACITY 50 lt

* External case in steel painted with epoxy paint ALL CAPACITY

* Stainless steel tank (made by mould)

* General switch

* Digital thermoregulator

* Safety thermostat

* Good insulation between external case and tank



Features

Item	610.1006.24
Temperature	+5 amb. to +120
Precision	°C ±0.5
Heating power	W 2100
Internal dimension (LxWxH)	mm 255x440x150
External dimension (LxWxH)	mm 445x560x240
Power	V/Hz 230/50-60
Weight	Kg 9
Security class	3.1
Protection class	IP 42
Capacity	lt 18.5
Unload + level	Si

