



Digital water bath with rectangular tank 50 lt - WB-MD50

Summary

*CAPACITY 5 - 15 - 24 lt * External case in ABS CAPACITY 50 lt * External case in steel painted with epoxy paint ALL CAPACITY

Description

*CAPACITY 5 - 15 - 24 lt

* External case in ABS CAPACITY 50 lt

* External case in steel painted with epoxy paint ALL CAPACITY

* Stainless steel tank (made by mould)

* General switch

* Digital thermoregulator

* Safety thermostat

* Good insulation between external case and tank



Features

Item	610.1006.50
Temperature	+5 amb. to +120
Precision	°C ±1
Heating power	W 3000
Internal dimension (LxWxH)	mm 440x560x200
External dimension (LxWxH)	mm 835x670x320
Power	V/Hz 230/50-60
Weight	Kg 33
Security class	3.1
Protection class	IP 42
Capacity	lt 50
Unload + level	Si

