



---

## Digital water bath with squared tank 5 lt - WB-MD5

### Summary

\*CAPACITY 5 - 15 - 24 lt \* External case in ABS CAPACITY 50 lt \* External case in steel painted with epoxy paint ALL CAPACITY

### Description

\*CAPACITY 5 - 15 - 24 lt

\* External case in ABS CAPACITY 50 lt

\* External case in steel painted with epoxy paint ALL CAPACITY

\* Stainless steel tank (made by mould)

\* General switch

\* Digital thermoregulator

\* Safety thermostat

\* Good insulation between external case and tank



## Features

Item	610.1006.05
Temperature	+5 amb. to +120
Precision	°C ±0.5
Heating power	W 700
Internal dimension (LxWxH)	mm 100x250x150
External dimension (LxWxH)	mm 280x360x240
Power	V/Hz 230/50-60
Weight	Kg 5.2
Security class	3.1
Protection class	IP 42
Capacity	lt 5
Unload + level	No

